How to Use the Handihaler



Instructions:

- 1. Open the cap by pulling upward.
- 2. Open the mouthpiece.
- 3. Peel back the foil from the medication capsule, remove the capsule and place the capsule in the Handihaler inhaler chamber.
- 4. Close the mouthpiece until you hear a click.
- 5. Press the green piercing button in once and release. This will release the medication from the capsule into the inhaler chamber. Do not shake the inhaler.
- 6. Turn your head away from the inhaler and exhale.
- 7. Place the mouthpiece between your teeth without biting it. Close your lips around the mouthpiece to form a good seal. Breathe in slowly and deeply so the capsule vibrates.
- 8. Remove the inhaler from your mouth.
- 9. Hold your breath for 10 seconds or as long as comfortable.
- 10. Exhale slowly through your mouth or nose away from the inhaler.
- 11. Repeat steps 7-9 to take the full dose.
- 12. Open mouthpiece and remove the used capsule.
- 13. Close the mouthpiece.
- 14. Close the cap.

Tips:

- If you are using a corticosteroid medication, rinse your mouth after all doses are complete.
- Always close the inhaler after use.