

CADD PRIZM-Intermittent

Equipment

- ✓ Medication at room temperature with attached tubing
- ✓ Pump
- ✓ 9 volt battery
- ✓ Saline flush
- ✓ Alcohol wipes
- ✓ Coin
- ✓ Sharps container
- ✓ Waste container

Procedure

1. **Wash your Hands**
2. **Clean your Work Area**
3. **Gather up your Equipment**
4. **Check your IV site (check for signs of redness, tenderness, swelling or drainage-notify your nurse if any signs present)**
5. **Check your Medication for:**
 - The right name
 - The correct medication and dose: _____
 - Particles or discoloration
 - Expiration date
6. **Stop the Pump**
 - Press the **“START/STOP”** pump
 - **“Stop the pump?”** will appear on the screen
 - Press the **“Y”** button
7. **Disconnect Completed Medication**
 - Unscrew tubing from your IV
8. **Administer Saline Flush**
 - Clean the cap on your IV with an alcohol wipe
 - Remove the cap from the saline flush
 - Remove any air from syringe
 - Twist the saline syringe onto your IV
 - Depress the plunger slowly to administer the flush
 - Remove the syringe and discard

9. Remove Medication/Cassette

- Close the clamp on the used medication tubing
- Unlatch the medication from the pump using a coin to release the lock
- Remove and discard the used medication

10. Change the Battery

- Slide the battery cover on the side of the pump
- Remove the old battery
- Insert new battery
- Slide cover back into the pump until closed
- Set pump down and wait for pump to complete self test

11. Attach new Medication Cassette or Bag

- Insert the new medication reservoir hooks into the hinge pins on the pump
- Place the pump upright on a flat surface
- Using the coin, push in and turn the latching button until the mark on the latch lines up with the solid dot and you hear 1 beep

12. Reset the Pump Program

- The pump will ask **“Reset Reservoir Volume to _____ ml?”**
- Press the **“Y”** button
- The pump will ask **“Prime tubing?”**
- Press the **“Y”** or **“N”** button
- Press and hold **“Y”** button until a few drops come out of end of tubing if you do need to prime

13. Reconnect Medication

- Clean the end of your IV with an alcohol wipe
- Remove the cap from the end of the tubing
- Twist the tubing onto your IV
- Open all clamps

14. Start the Pump

- Press the **“START/STOP”** button
- The pump will ask **“Start the Pump?”**
- Press the **“Y”** button

The pump will run through the program. The Green Light will begin to flash.

Your Medication is due every _____ hours.

The times your pump will administer your medication are _____.

Change your medication bag/cassette every _____ hours, between the hours of _____ and _____.